

Explore our newly expanded range of Grignard Reagents

The Grignard reaction is an essential method in the formation of carbon-carbon bonds.

In response to customer demand for more efficient Grignard Reagents and the desire to carry out chemistry in an environmentally friendly manner, we recently expanded our range of Grignard Reagent products by almost 50%. Supplied in solvents such as diethyl ether, tetrahydrofuran, 2-methyltetrahydrofuran and cyclopropylmethyl ether, our Grignard Reagent range offers the widest choice of reactivity, ease of work up and reduced environmental impact.

In reactions involving Grignard Reagents it is essential to exclude water and air, which rapidly destroy the reagent. Therefore, we provide these sensitive reagents in the industry-leading AcroSeal® packaging, preventing the degradation of the reagents.

We offer the following Grignard Reagents.

* New  Green Solvent

Product Code	Product Name	Packaging Size	CASNo
209531000	Allylmagnesium bromide, 1M solution in diethyl ether, AcroSeal	100mL	1730-25-2
209538000	Allylmagnesium bromide, 1M solution in diethyl ether, AcroSeal	800mL	1730-25-2
209671000	Allylmagnesium chloride, 1.7M solution in THF, AcroSeal	100mL	2622-05-1
209678000	Allylmagnesium chloride, 1.7M solution in THF, AcroSeal	800mL	2622-05-1
331651000	Benzylmagnesium chloride, 1.4M solution in THF, AcroSeal	100mL	6921-34-2
331658000	Benzylmagnesium chloride, 1.4M solution in THF, AcroSeal	800mL	6921-34-2
389551000	Benzylmagnesium chloride, 1.5M solution in THF, AcroSeal	100mL	6921-34-2
389558000	Benzylmagnesium chloride, 1.5M solution in THF, AcroSeal	800mL	6921-34-2
377391000	Benzylmagnesium chloride, 1M solution in diethyl ether, AcroSeal	100mL	6921-34-2
377398000	Benzylmagnesium chloride, 1M solution in diethyl ether, AcroSeal	800mL	6921-34-2
317371000	Benzylmagnesium chloride, 2M solution in THF, AcroSeal	100mL	6921-34-2
317378000	Benzylmagnesium chloride, 2M solution in THF, AcroSeal	800mL	6921-34-2
440320500	2-Benzyloxyphenylmagnesium bromide, 1M solution in THF, AcroSeal	50mL	328000-16-4
440330500	3-Benzyloxyphenylmagnesium bromide, 1M solution in THF, AcroSeal	50mL	36281-96-6
436221000	4-Benzyloxyphenylmagnesium bromide, 0.8M solution in THF, AcroSeal	100mL	120186-59-6
431621000	4-Biphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	3315-91-1
353250010	Bis(cyclopentadienyl)magnesium, sublimed	1g	1284-72-6
353252500	Bis(cyclopentadienyl)magnesium, sublimed	250mg	1284-72-6
392151000	3,5-Bis(trifluoromethyl)phenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	112981-69-8
431611000	2-Bromobenzylmagnesium bromide, 0.25M solution in diethyl ether, AcroSeal	100mL	56812-60-3
431601000	3-Bromobenzylmagnesium bromide, 0.25M solution in diethyl ether, AcroSeal	100mL	107549-22-4
427461000	3-Butenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	7103-09-5
427468000	3-Butenylmagnesium bromide, 0.5M solution in THF, AcroSeal	800mL	7103-09-5
438711000	2-Butenylmagnesium chloride, 0.5M solution in THF, AcroSeal	100mL	22649-70-3
438718000	2-Butenylmagnesium chloride, 0.5M solution in THF, AcroSeal	800mL	22649-70-3
381181000	n-Butylethylmagnesium, 0.9M solution in heptane, AcroSeal	100mL	62202-86-2
381188000	n-Butylethylmagnesium, 0.9M solution in heptane, AcroSeal	800mL	62202-86-2
377781000	n-Butylethylmagnesium, 1.3M solution in heptane, AcroSeal	100mL	62202-86-2
377788000	n-Butylethylmagnesium, 1.3M solution in heptane, AcroSeal	800mL	62202-86-2
394301000	tert-Butylmagnesium chloride, 1.7M solution in diethyl ether, AcroSeal	100mL	677-22-5
394308000	tert-Butylmagnesium chloride, 1.7M solution in diethyl ether, AcroSeal	800mL	677-22-5
331671000	tert-Butylmagnesium chloride, 1.7M solution in THF, AcroSeal	100mL	677-22-5
331678000	tert-Butylmagnesium chloride, 1.7M solution in THF, AcroSeal	800mL	677-22-5
208541000	tert-Butylmagnesium chloride, 2.0M solution in diethyl ether, AcroSeal	100mL	677-22-5
208548000	tert-Butylmagnesium chloride, 2.0M solution in diethyl ether, AcroSeal	800mL	677-22-5
331661000	sec-Butylmagnesium chloride, 25 wt% solution in THF, AcroSeal	100mL	15366-08-2
331681000	n-Butylmagnesium chloride, 20 wt% solution in THF/toluene, AcroSeal	100mL	693-04-9
331688000	n-Butylmagnesium chloride, 20 wt% solution in THF/toluene, AcroSeal	800mL	693-04-9
432861000	4-tert-Butylphenylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	100mL	63488-10-8
426971000	n-Butyl-sec-butylmagnesium, 0.7M solution in hexane, AcroSeal	100mL	39881-32-8
431990500	3-Chloro-4-fluorophenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	413589-34-1
432005000	3-Chloro-5-fluorophenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	480438-50-4

Product Code	Product Name	Packaging Size	CASNo
399081000	4-Chloro-3-fluorophenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	170793-00-7
431980500	4-Chloro-2-methylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	480438-47-9
397490500	3-Chlorophenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	36229-42-2
397491000	3-Chlorophenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	36229-42-2
377401000	4-Chlorophenylmagnesium bromide, 1M solution in diethyl ether, AcroSeal	100mL	873-77-8
377408000	4-Chlorophenylmagnesium bromide, 1M solution in diethyl ether, AcroSeal	800mL	873-77-8
377251000	4-Chlorophenylmagnesium bromide, 1M solution in THF/toluene, AcroSeal	100mL	873-77-8
377258000	4-Chlorophenylmagnesium bromide, 1M solution in THF/toluene, AcroSeal	800mL	873-77-8
431970500	5-Chloro-2-thienylmagnesium bromide, 0.5M slurry in THF, AcroSeal	50mL	111762-30-2
377341000	Cyclohexylmagnesium chloride, 1.3M solution in THF/toluene, AcroSeal	100mL	931-51-1
377348000	Cyclohexylmagnesium chloride, 1.3M solution in THF/toluene, AcroSeal	800mL	931-51-1
433611000	Cyclohexylmethylmagnesium bromide, 0.25M solution in THF, AcroSeal	100mL	35166-78-0
427451000	Cyclopentylmagnesium bromide, 2.0M solution in diethyl ether, AcroSeal	100mL	33240-34-5
427458000	Cyclopentylmagnesium bromide, 2.0M solution in diethyl ether, AcroSeal	800mL	33240-34-5
208461000	Cyclopentylmagnesium chloride, 2M solution in diethyl ether, AcroSeal	100mL	32916-51-1
208468000	Cyclopentylmagnesium chloride, 2M solution in diethyl ether, AcroSeal	800mL	32916-51-1
397611000	Cyclopropylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	23719-80-4
431570500	3,4-Dichlorophenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	79175-35-2
433010500	3,5-Dichlorophenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	82297-90-3
399511000	3,4-Difluorophenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	90897-92-0
427441000	3,5-Difluorophenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	62351-47-7
431551000	3,5-Dimethoxybenzylmagnesium chloride, 0.2M solution in THF, AcroSeal	100mL	135808-66-1
431960500	3,5-Dimethyl-4-methoxyphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	185416-17-5
431930500	2,3-Dimethylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	134640-85-0
431950500	2,4-Dimethylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	34589-46-3
431920500	2,5-Dimethylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	30897-86-0
436061000	2,6-Dimethylphenylmagnesium bromide, 1M solution in THF, AcroSeal	100mL	21450-64-6
436068000	2,6-Dimethylphenylmagnesium bromide, 1M solution in THF, AcroSeal	800mL	21450-64-6
431910500	3,4-Dimethylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	89980-68-7
431940500	3,5-Dimethylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	34696-73-6
436051000	1,1-Dimethylpropylmagnesium chloride, 1M solution in diethyl ether, AcroSeal	100mL	28276-08-6
427431000	2,2-Dimethylpropylmagnesium chloride, 1.0M solution in diethyl ether, AcroSeal	100mL	13132-23-5
377771000	Di-n-butylmagnesium, 0.5M solution in heptane, AcroSeal	100mL	1191-47-5
377778000	Di-n-butylmagnesium, 0.5M solution in heptane, AcroSeal	800mL	1191-47-5
431900500	(1,3-Dioxan-2-ylethyl)magnesium bromide, 0.5M solution in THF, AcroSeal	50mL	480438-44-6
438861000	(1,3-Dioxolan-2-ylmethyl)magnesium bromide, 0.5M solution in THF, AcroSeal	100mL	180675-22-3
438868000	(1,3-Dioxolan-2-ylmethyl)magnesium bromide, 0.5M solution in THF, AcroSeal	800mL	180675-22-3
435951000	(2-Ethylhexyl)magnesium bromide, 1M solution in diethyl ether, AcroSeal	100mL	90224-21-8
347281000	Ethylmagnesium bromide, 0.9M in tert-Butyl methyl ether, AcroSeal	100mL	925-90-6
347288000	Ethylmagnesium bromide, 0.9M in tert-Butyl methyl ether, AcroSeal	800mL	925-90-6
210471000	Ethylmagnesium bromide, 0.9M solution in THF, AcroSeal	100mL	925-90-6
210478000	Ethylmagnesium bromide, 0.9M solution in THF, AcroSeal	800mL	925-90-6
347290010	Ethylmagnesium bromide, 3M in diethyl ether, AcroSeal	1L	925-90-6
347291000	Ethylmagnesium bromide, 3M in diethyl ether, AcroSeal	100mL	925-90-6
347298000	Ethylmagnesium bromide, 3M in diethyl ether, AcroSeal	800mL	925-90-6
183531000	Ethylmagnesium bromide, ca. 3.2M solution in 2-MeTHF; AcroSeal	100mL	925-90-6
183538000	Ethylmagnesium bromide, ca. 3.2M solution in 2-MeTHF; AcroSeal	800mL	925-90-6
377281000	Ethylmagnesium chloride, 2.0M solution in THF, AcroSeal	100mL	2386-64-3
377288000	Ethylmagnesium chloride, 2.0M solution in THF, AcroSeal	800mL	2386-64-3
252570024	Ethylmagnesium chloride, 2.7M (25 wt.%) solution in THF, AcroSeal	2.4L	2386-64-3
252571000	Ethylmagnesium chloride, 2.7M (25 wt.%) solution in THF, AcroSeal	100mL	2386-64-3
252578000	Ethylmagnesium chloride, 2.7M (25 wt.%) solution in THF, AcroSeal	800mL	2386-64-3
431890500	4-Ethylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	22873-28-5
439121000	Ethyneylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	4301-14-8
439128000	Ethyneylmagnesium bromide, 0.5M solution in THF, AcroSeal	800mL	4301-14-8
388951000	Ethyneylmagnesium chloride, 0.6M solution in THF/toluene, AcroSeal	100mL	65032-27-1
388958000	Ethyneylmagnesium chloride, 0.6M solution in THF/toluene, AcroSeal	800mL	65032-27-1
377441000	Ethyneylmagnesium chloride, 1.1M solution in THF/toluene, AcroSeal	100mL	65032-27-1
377448000	Ethyneylmagnesium chloride, 1.1M solution in THF/toluene, AcroSeal	800mL	65032-27-1
426790500	4-Fluorobenzylmagnesium chloride, 0.25M solution in THF, AcroSeal	50mL	1643-73-8
431870500	3-Fluoro-4-methoxyphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	112780-67-3
431880500	5-Fluoro-2-methoxyphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	188132-02-7
433740500	3-Fluoro-2-methylphenylmagnesium chloride, 0.5M solution in THF, AcroSeal	50mL	480424-74-6
431850500	3-Fluoro-4-methylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	185077-02-5
431840500	4-Fluoro-2-methylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	30897-90-6
431860500	4-Fluoro-3-methylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	82297-89-0
431830500	5-Fluoro-2-methylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	186496-59-3
427821000	3-Fluorophenylmagnesium bromide, 1M solution in THF, AcroSeal	100mL	17318-03-5
377231000	4-Fluorophenylmagnesium bromide, 0.8M solution in THF, AcroSeal	100mL	352-13-6
377238000	4-Fluorophenylmagnesium bromide, 0.8M solution in THF, AcroSeal	800mL	352-13-6
377371000	4-Fluorophenylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	100mL	352-13-6
377378000	4-Fluorophenylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	800mL	352-13-6
433031000	Heptylmagnesium bromide, 1M solution in diethyl ether, AcroSeal	100mL	13125-66-1
427761000	Hexylmagnesium bromide, 0.8M solution in THF, AcroSeal	100mL	3761-92-0
427768000	Hexylmagnesium bromide, 0.8M solution in THF, AcroSeal	800mL	3761-92-0
434681000	Hexylmagnesium chloride, 2M solution in THF, AcroSeal	100mL	44767-62-6
434688000	Hexylmagnesium chloride, 2M solution in THF, AcroSeal	800mL	44767-62-6
427391000	Isobutylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	100mL	926-62-5
427398000	Isobutylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	800mL	926-62-5

Product Code	Product Name	Packaging Size	CASNo
209521000	Isobutylmagnesium chloride, 2M solution in diethyl ether, AcroSeal	100mL	5674-02-2
209528000	Isobutylmagnesium chloride, 2M solution in diethyl ether, AcroSeal	800mL	5674-02-2
427751000	Isopropenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	13291-18-4
427758000	Isopropenylmagnesium bromide, 0.5M solution in THF, AcroSeal	800mL	13291-18-4
426781000	Isopropylmagnesium bromide, 3M solution in 2-MeTHF, AcroSeal	100mL	920-39-8
386281000	Isopropylmagnesium chloride - lithium chloride complex, 1.3M solution in THF, AcroSeal	100mL	745038-86-2
386288000	Isopropylmagnesium chloride - lithium chloride complex, 1.3M solution in THF, AcroSeal	800mL	745038-86-2
212851000	Isopropylmagnesium chloride, 2.0M solution in THF, AcroSeal	100mL	1068-55-9
212858000	Isopropylmagnesium chloride, 2.0M solution in THF, AcroSeal	800mL	1068-55-9
431630500	4-Isopropylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	18620-03-6
413450250	Magnesium methyl carbonate, pure, 1.8M solution in DMF	25g	4861-79-4
210731000	2-Mesitylmagnesium bromide, 1M solution in THF, AcroSeal	100mL	2633-66-1
433730500	2-Methoxybenzylmagnesium chloride, 0.25M solution in THF, AcroSeal	50mL	480438-46-8
431820500	3-Methoxybenzylmagnesium chloride, 0.25M solution in THF, AcroSeal	50mL	26905-40-8
427421000	4-Methoxybenzylmagnesium chloride, 0.25M solution in THF, AcroSeal	100mL	38769-92-5
431810500	4-Methoxy-2-methylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	148706-30-3
431800500	6-Methoxy-2-naphthylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	38046-82-1
440800500	4-Methoxyphenethylmagnesium chloride, 0.5M solution in THF, AcroSeal	50mL	211115-05-8
426730010	2-Methoxyphenylmagnesium bromide, 1M solution in THF, AcroSeal	1L	16750-63-3
426731000	2-Methoxyphenylmagnesium bromide, 1M solution in THF, AcroSeal	100mL	16750-63-3
426738000	2-Methoxyphenylmagnesium bromide, 1M solution in THF, AcroSeal	800mL	16750-63-3
426751000	3-Methoxyphenylmagnesium bromide, 1M solution in THF/toluene, AcroSeal	100mL	36282-40-3
426758000	3-Methoxyphenylmagnesium bromide, 1M solution in THF/toluene, AcroSeal	800mL	36282-40-3
377421000	4-Methoxyphenylmagnesium bromide, 1M solution in THF, AcroSeal	100mL	13139-86-1
377428000	4-Methoxyphenylmagnesium bromide, 1M solution in THF, AcroSeal	800mL	13139-86-1
427411000	2-Methylallylmagnesium chloride, 0.5M solution in THF, AcroSeal	100mL	5674-01-1
433710500	2-Methylbenzylmagnesium chloride, 0.25M solution in THF, AcroSeal	50mL	29875-05-6
433720500	3-Methylbenzylmagnesium chloride, 0.25M solution in THF, AcroSeal	50mL	29875-06-7
433700500	4-Methylbenzylmagnesium chloride, 0.5M solution in THF, AcroSeal	50mL	29875-07-8
439761000	3,4-(Methylenedioxy)phenylmagnesium bromide, 0.8M solution in toluene/THF (50:50), AcroSeal	100mL	17680-04-5
377331000	Methylmagnesium bromide, 1.4M solution in toluene/THF, AcroSeal	100mL	75-16-1
377338000	Methylmagnesium bromide, 1.4M solution in toluene/THF, AcroSeal	800mL	75-16-1
377381000	Methylmagnesium bromide, 1M solution in THF, AcroSeal	100mL	75-16-1
377388000	Methylmagnesium bromide, 1M solution in THF, AcroSeal	800mL	75-16-1
391121000	Methylmagnesium bromide, 3.2M solution in 2-MeTHF, AcroSeal	100mL	75-16-1
391128000	Methylmagnesium bromide, 3.2M solution in 2-MeTHF, AcroSeal	800mL	75-16-1
183541000	Methylmagnesium bromide, 3M solution in diethyl ether, AcroSeal	100mL	75-16-1
183548000	Methylmagnesium bromide, 3M solution in diethyl ether, AcroSeal	800mL	75-16-1
252561000	Methylmagnesium chloride, 3M (22 wt.%) solution in THF, AcroSeal	100mL	676-58-4
252568000	Methylmagnesium chloride, 3M (22 wt.%) solution in THF, AcroSeal	800mL	676-58-4
427401000	Methylmagnesium iodide, 3M solution in diethyl ether, AcroSeal	100mL	917-64-6
427408000	Methylmagnesium iodide, 3M solution in diethyl ether, AcroSeal	800mL	917-64-6
433680500	2-Methyl-1-naphthylmagnesium bromide, 0.25M solution in THF, AcroSeal	50mL	21450-90-8
433690500	4-Methyl-1-naphthylmagnesium bromide, 0.25M solution in THF, AcroSeal	50mL	115913-58-1
438751000	2-Methyl-2-phenylpropylmagnesium chloride, 0.5M solution in diethyl ether, AcroSeal	100mL	35293-35-7
438758000	2-Methyl-2-phenylpropylmagnesium chloride, 0.5M solution in diethyl ether, AcroSeal	800mL	35293-35-7
436021000	1-Methyl-1-propenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	85676-85-3
436028000	1-Methyl-1-propenylmagnesium bromide, 0.5M solution in THF, AcroSeal	800mL	85676-85-3
435561000	2-Methyl-1-propenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	38614-36-7
435568000	2-Methyl-1-propenylmagnesium bromide, 0.5M solution in THF, AcroSeal	800mL	38614-36-7
431780500	3-Methyl-2-thienylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	95184-07-9
431750500	[2-(4-Morpholinylmethyl)phenyl]magnesium bromide, 0.25M solution in THF, AcroSeal	50mL	480424-77-9
431760500	[3-(4-Morpholinylmethyl)phenyl]magnesium bromide, 0.25M solution in THF, AcroSeal	50mL	480424-76-8
431770500	[4-(4-Morpholinylmethyl)phenyl]magnesium bromide, 0.25M solution in THF, AcroSeal	50mL	480424-75-7
440791000	(2-Naphthalenylmethyl)magnesium bromide, 0.25M solution in diethyl ether, AcroSeal	100mL	127543-80-0
384120500	1-Naphthylmagnesium bromide, 0.25M solution in THF, AcroSeal	50mL	703-55-9
384121000	1-Naphthylmagnesium bromide, 0.25M solution in THF, AcroSeal	100mL	703-55-9
431740500	2-Naphthylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	21473-01-8
440781000	Nonylmagnesium bromide, 1M solution in diethyl ether, AcroSeal	100mL	39691-62-8
436011000	Octadecylmagnesium chloride, 0.3M solution in THF, AcroSeal	100mL	116980-66-6
436018000	Octadecylmagnesium chloride, 0.3M solution in THF, AcroSeal	800mL	116980-66-6
427811000	Octylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	100mL	17049-49-9
427818000	Octylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	800mL	17049-49-9
428031000	Octylmagnesium chloride, 1.4M solution in THF, AcroSeal	100mL	38841-98-4
428038000	Octylmagnesium chloride, 1.4M solution in THF, AcroSeal	800mL	38841-98-4
436160500	Pentadecylmagnesium bromide, 0.5M solution in diethyl ether, AcroSeal	50mL	78887-70-4
438741000	Pentafluorophenylmagnesium bromide, 0.5M solution in diethyl ether, AcroSeal	100mL	879-05-0
436001000	Pentamethylenebis(magnesium bromide), 0.5M solution in THF, AcroSeal	100mL	23708-48-7
435961000	3-Pentylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	100mL	4852-26-0
435551000	Pentylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	100mL	693-25-4
435558000	Pentylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	800mL	693-25-4
436211000	Pentylmagnesium chloride, 1.3M solution in THF, AcroSeal	100mL	6393-56-2
436218000	Pentylmagnesium chloride, 1.3M solution in THF, AcroSeal	800mL	6393-56-2
431730500	9-Phenanthrylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	71112-64-6
427801000	Phenethylmagnesium chloride, 1M solution in THF, AcroSeal	100mL	90878-19-6
427808000	Phenethylmagnesium chloride, 1M solution in THF, AcroSeal	800mL	90878-19-6
431720500	4-Phenoxyphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	21473-02-9
377411000	Phenylmagnesium bromide, 1M solution in THF, AcroSeal	100mL	100-58-3
377418000	Phenylmagnesium bromide, 1M solution in THF, AcroSeal	800mL	100-58-3

Product Code	Product Name	Packaging Size	CASNo
391111000	Phenylmagnesium bromide, 2.8M solution in 2-MeTHF, AcroSeal	100mL	100-58-3
391118000	Phenylmagnesium bromide, 2.8M solution in 2-MeTHF, AcroSeal	800mL	100-58-3
173981000	Phenylmagnesium bromide, 2.8M solution in diethyl ether, AcroSeal	100mL	100-58-3
173988000	Phenylmagnesium bromide, 2.8M solution in diethyl ether, AcroSeal	800mL	100-58-3
252581000	Phenylmagnesium chloride, 25 wt% solution in THF, AcroSeal	100mL	100-59-4
252588000	Phenylmagnesium chloride, 25 wt% solution in THF, AcroSeal	800mL	100-59-4
431710500	[2-(1-Piperidinylmethyl)phenyl]magnesium bromide, 0.25M solution in THF, AcroSeal	50mL	480424-81-5
431700500	[3-(1-Piperidinylmethyl)phenyl]magnesium bromide, 0.25M solution in THF, AcroSeal	50mL	480424-82-6
431690500	[4-(1-Piperidinylmethyl)phenyl]magnesium bromide, 0.25M solution in THF, AcroSeal	50mL	480424-83-7
434671000	1-Propenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	14092-04-7
434678000	1-Propenylmagnesium bromide, 0.5M solution in THF, AcroSeal	800mL	14092-04-7
208451000	Propylmagnesium chloride, 2M solution in diethyl ether, AcroSeal	100mL	2234-82-4
208458000	Propylmagnesium chloride, 2M solution in diethyl ether, AcroSeal	800mL	2234-82-4
431680500	4-n-Propylphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	87942-08-3
426071000	1-Propynylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	16466-97-0
426078000	1-Propynylmagnesium bromide, 0.5M solution in THF, AcroSeal	800mL	16466-97-0
431650500	[2-(1-Pyrrolidinylmethyl)phenyl]magnesium bromide, 0.25M solution in THF, AcroSeal	50mL	480424-80-4
431660500	[3-(1-Pyrrolidinylmethyl)phenyl]magnesium bromide, 0.25M solution in THF, AcroSeal	50mL	480424-79-1
431670500	[4-(1-Pyrrolidinylmethyl)phenyl]magnesium bromide, 0.25M solution in THF, AcroSeal	50mL	480424-78-0
431640500	4-(2-Tetrahydro-2H-pyranoyl)phenylmagnesium bromide, 0.5M solution in THF, AcroSeal	50mL	36637-44-2
434610500	2-Thienylmagnesium bromide, 1M solution in THF, AcroSeal	50mL	5713-61-1
433650500	3-Thienylmagnesium iodide, 0.3M solution in THF, AcroSeal	50mL	172035-86-8
399071000	4-Thioanisolemagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	18620-04-7
427471000	o-Tolylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	100mL	932-31-0
427478000	o-Tolylmagnesium bromide, 2M solution in diethyl ether, AcroSeal	800mL	932-31-0
426761000	p-Tolylmagnesium bromide, ca. 0.4M(10%) solution in diethyl ether, AcroSeal	100mL	4294-57-9
426768000	p-Tolylmagnesium bromide, ca. 0.4M(10%) solution in diethyl ether, AcroSeal	800mL	4294-57-9
435971000	m-Tolylmagnesium chloride, 1M solution in THF, AcroSeal	100mL	121905-60-0
435978000	m-Tolylmagnesium chloride, 1M solution in THF, AcroSeal	800mL	121905-60-0
427481000	o-Tolylmagnesium chloride, 1.4M (22 wt.%) solution in THF/toluene, AcroSeal	100mL	33872-80-9
427488000	o-Tolylmagnesium chloride, 1.4M (22 wt.%) solution in THF/toluene, AcroSeal	800mL	33872-80-9
382721000	p-Tolylmagnesium chloride, 2M (31 wt.%) solution in THF, AcroSeal	100mL	696-61-7
382728000	p-Tolylmagnesium chloride, 2M (31 wt.%) solution in THF, AcroSeal	800mL	696-61-7
433601000	3,4,5-Trifluorophenylmagnesium bromide, 0.3M solution in THF, AcroSeal	100mL	156006-28-9
431541000	3,4,5-Trimethoxyphenylmagnesium bromide, 0.5M solution in THF, AcroSeal	100mL	133095-91-7
377461000	(Trimethylsilyl)methylmagnesium chloride, 1.3M solution in THF, AcroSeal	100mL	13170-43-9
377468000	(Trimethylsilyl)methylmagnesium chloride, 1.3M solution in THF, AcroSeal	800mL	13170-43-9
209391000	Vinylimagnesium bromide, 0.7M solution in THF, AcroSeal	100mL	1826-67-1
209398000	Vinylimagnesium bromide, 0.7M solution in THF, AcroSeal	800mL	1826-67-1
252591000	Vinylimagnesium chloride, 1.9M (16.5 wt.%) solution in THF, AcroSeal	100mL	3536-96-7
252598000	Vinylimagnesium chloride, 1.9M (16.5 wt.%) solution in THF, AcroSeal	800mL	3536-96-7



Products are available in a variety of pack sizes. For more information please visit our web site www.acros.com or contact your local distributor.

Austria	Belgium	Czech Republic	Denmark	Finland	France
Fisher Scientific (Austria) GmbH Rudolf Alt-Platz 1 A-1030 Wien Toll-Free Number: 0800 20 88 40 Fax: 0800 20 66 90 info.austria@thermofisher.com www.at.fishersci.com	Fisher Scientific Clintonpark - Keppekouter Ninovesteenweg 198 B-9320 Erembodegem Tel: +32 (0) 56 260 260 Fax: +32 (0) 56 260 270 be.fisher@thermofisher.com www.be.fishersci.com	Fisher Scientific, spol. s r.o. Kosmonauta 324 CZ-530 09 Pardubice Tel: +420 466 798 230 Fax: +420 466 435 008 info.cz@thermofisher.com www.thermofisher.cz	Fisher Scientific Industrivej 3 DK-3550 Slangerup Tel: +45 70 27 99 20 Fax: +45 70 27 99 29 kundeservice@thermofisher.com www.fishersci.dk	Fisher Scientific Oy Ratastie 2 FI-01620 Vantaa Tel: +358 802 76 280 Fax: +358 802 76 235 fisher.fi@thermofisher.com www.fishersci.fi	Fisher Scientific FAS Parc d'Innovation BP 50111 67403 Illkirch Cedex Tel: +33 (0) 3 88 67 53 20 Fax: +33 (0) 3 88 67 11 68 fr.commande@thermofisher.com www.fr.fishersci.com
Germany	Ireland	Italy	The Netherlands	Norway	Portugal
Fisher Scientific GmbH Im Heiligen Feld 17 D-58239 Schwerte Toll-Free Number: 0800 3 47 43 70 Fax: 0800 3 47 43 71 info.germany@thermofisher.com www.de.fishersci.com	Fisher Scientific Ireland Suite 3 Plaza 212 Blanchardstown Corporate Park 2 Ballycoolin Dublin 15 Tel: +353 01 885 5854 Fax: +353 01 899 1855 fsi.sales@thermofisher.com www.ie.fishersci.com	Fisher Scientific Tel: +39 02 953 28 258 Fax: +39 02 953 27 374 it.fisher@thermofisher.com www.it.fishersci.com	Fisher Scientific Postbus 4 Scheepsbouwersweg 1b 1120 AA Landsmeer Tel: +31 (0) 20 487 70 70 Fax: +31 (0) 20 487 70 70 nl.info@thermofisher.com www.fishersci.nl	Fisher Scientific Frysjaeven 33E 0884 Oslo Tel: +47 22 95 59 59 Fax: +47 22 95 59 40 fisher.no@thermofisher.com www.fishersci.no	Fisher Scientific Rua Pedro Álvares Cabral, n°24, 3ºD Edifício Euro - Infantado 2670-391 Loures Tel: +351 21 425 33 50/4 Fax: +351 21 425 33 51 pt.fisher@thermofisher.com www.pt.fishersci.com
Slovenia	Spain	Sweden	Switzerland	United Kingdom	
Fisher Scientific d.o.o. Parmova 53 SI-1000 Ljubljana Tel: +386 597 31 494 Fax: +386 597 31 495 si.info@thermofisher.com www.si.fishersci.com	Fisher Scientific C/Luis I, 9 28031 Madrid Tel: +34 91 380 67 10 Fax: +34 91 380 85 02 es.fisher@thermofisher.com www.es.fishersci.com	Fisher Scientific Box 9193 400 94 Göteborg Tel: +46 31 352 32 00 Fax: +46 31 352 32 50 fisher.se@thermofisher.com www.fishersci.se	Fisher Scientific AG Wilstrasse 57 CH-5610 Wohlen Tel: +41 (0) 56 618 41 11 Fax: +41 (0) 56 618 41 41 info.ch@thermofisher.com www.ch.fishersci.com	Fisher Scientific UK Ltd. Bishop Meadow Road Loughborough Leicestershire LE11 5RG Tel: +44 (0)1509 555500 Fax: +44 (0)1509 555111 fsuk.sales@thermofisher.com www.fisher.co.uk	

